Válido para todos os países, exceto Estados Unidos / Valid for all countries, except United States / Válido para los países, excepto Estados Unidos / Valido per tutti paesi, ad eccezione del Stati Uniti / Gültig für alle Länder, außer Vereinigte Staaten / Valable pour tous les pays, sauf le Etats-Unis / Gælder for alle lande, undtagen USA / Az Egyesült Államok kivételével minden országra érvényes / Valabil pentru toate țările, cu excepția Statelor Unite / Giltig i alla länder utanför USA

Válido apenas para EUA / Valid only for the United States / Válido solamente para Estados Unidos / Valido solo per gli Stati Uniti / Gültig nur für US / Valide seulement pour les Etats-Unis / Gælder kun for USA / Kizárólag az Egyesült Államokra érvényes / Valabil numai pentru Statele Unite / Giltig endast i USA

Português	Base de Titânio
English	Titanium Base
Español	Base de Titanio
Italiano	Base in Titanio
Deutsch	Titanbasis
Français	Base en Titane
Dansk	Titanbase
Magyar	Titán bázis
Română	Bază de titan
Svenska	Titanbas

English	Titanium Base, US only
---------	------------------------

135.009; 135.009-ET; 135.010; 135.010-ET; 135.011; 135.011-ET; 135.012; 135.012-ET; 135.017; 135.017-ET; 135.018; 135.018-ET; 135.019; 135.019-ET; 135.020; 135.020-ET; 135.025; 135.025-ET; 135.026; 135.026-ET; 135.027; 135.027-ET; 135.028; 135.028-ET; 135.059; 135.060; 135.061; 135.062; 135.067; 135.068; 135.069; 135.070; 135.075; 135.076; 135.076-ET; 135.077; 135.077-ET; 135.078; 135.078-ET; 135.083; 135.084; 135.085; 135.086; 135.101; 135.102; 135.103; 135.104; 135.113; 135.113-ET; 135.114; 135.114-ET; 135.115; 135.115-ET; 135.116; 135.116-ET; 135.117; 135.117-ET; 135.202; 135.202-ET; 135.203; 135.203-ET; 135.204; 135.204-ET; 135.205; 135.205-ET; 135.206; 135.206-ET; 135.207; 135.208; 135.208-ET; 135.209; 135.209-ET; 135.210; 135.210-ET; 135.211; 135.211-ET; 135.212; 135.212-ET; 135.213; 135.214; 135.214-ET; 135.215; 135.215-ET; 135.216; 135.216-ET; 135.217; 135.217-ET; 135.218; 135.218-ET; 135.219; 135.220; 135.220-ET; 135.221; 135.221-ET; 135.222; 135.222-ET; 135.223; 135.223-ET; 135.224; 135.224-ET; 135.225; 135.225-ET; 135.236-ET; 135.237-ET; 135.238-ET; 135.239-ET; 135.241-ET; 135.242-ET; 135.243-ET; 135.244-ET; 135.245-ET; 135.260; 135.261; 135.262; 135.263; 135.264; 135.266; 135.267; 135.268; 135.269; 135.270; 135.272; 135.273; 135.274; 135.275; 135.276; 135.278; 135.279; 135.280; 135.281; 135.282; 135.284; 135.285; 135.286; 135.287; 135.288; 135.290; 135.291; 135.292; 135.293; 135.294; 135.304; 135.304-ET; 135.305; 135.305-ET; 135.306; 135.306-ET; 135.307; 135.307-ET; 135.308; 135.308-ET; 135.309; 135.309-ET; 135.310; 135.310-ET; 135.311; 135.311-ET; 135.312; 135.312-ET; 135.313; 135.313-ET; 135.314; 135.314-ET; 135.315; 135.315-ET; 135.316; 135.316-ET; 135.317; 135.317-ET; 135.318; 135.318-ET; 135.319; 135.319-ET; 135.320; 135.320-ET; 135.321; 135.321-ET; 135.322; 135.322-ET; 135.323; 135.323-ET; 135.324; 135.324-ET; 135.325; 135.325-ET; 135.326; 135.326-ET; 135.327; 135.327-ET; 135.328; 135.328-ET; 135.329; 135.329-ET; 135.330; 135.330-ET; 135.331; 135.331-ET; 135.332; 135.332-ET; 135.333; 135.333-ET; 135.334; 135.334-ET; 135.335; 135.335-ET; 135.336; 135.336-ET; 135.337; 135.337-ET; 135.338; 135.338-ET; 135.339; 135.339-ET; 135.340; 135.340-ET; 135.341; 135.341-ET; 135.342; 135.342-ET; 135.343; 135.343-ET; 135.344; 135.344-ET; 135.355; 135.355-ET; 135.356; 135.356-ET; 135.357; 135.357-ET; 135.358; 135.358-ET; 135.359; 135.359-ET; 135.361; 135.361-ET; 135.362; 135.362-ET; 135.363; 135.363-ET; 135.364; 135.364-ET; 135.365; 135.365-ET; 135.367; 135.367-ET; 135.368; 135.368-ET; 135.369; 135.369-ET; 135.370; 135.370-ET; 135.371; 135.371-ET; 135.373; 135.373-ET; 135.374; 135.374-ET; 135.375; 135.375-ET; 135.376; 135.376-ET; 135.377; 135.377-ET; 135.379; 135.379-ET; 135.380; 135.380-ET; 135.381; 135.381-ET; 135.382; 135.382-ET; 135.383; 135.383-ET; 135.385; 135.385-ET; 135.386; 135.386-ET; 135.387; 135.387-ET; 135.388; 135.388-ET; 135.389; 135.389-ET; 135.391; 135.391-ET; 135.392; 135.392-ET; 135.393; 135.393-ET; 135.394; 135.394-ET; 135.395; 135.395-ET; 135.396; 135.396-ET; 135.397; 135.397-ET; 135.398; 135.398-ET; 135.399; 135.399-ET; 135.400; 135.400-ET; 135.401; 135.401-ET; 135.402; 135.402-ET; 135.403; 135.403-ET; 135.404; 135.404-ET; 135.405; 135.405-ET; 135.406; 135.406-ET; 135.407; 135.407-ET; 135.408; 135.408-ET; 135.409; 135.409-ET; 135.410; 135.410-ET; 135.411; 135.411-ET; 135.412; 135.412-ET; 135.413; 135.413-ET

JJGC Indústria e Comércio de Materiais Dentários S.A. - Avenida Juscelino Kubitschek de Oliveira, 3.291 - Cidade Industrial - 81.270-200 - Curitiba - Paraná - Brasil - CNPJ : 00.489.050/0001-84
Nome fantasia / Brand name / Nombre de la marca / Nome commerciale / Markenname / Nom de marque / Mærkenavn / Márkanév / Numele mărcii / Varumärkesnamn : Neodent - Indústria Brasileira / Brazilian Industry / Industria Brasileira / Industria Brasileira / Brasilianische Industrie / Industrie brésilienne / Brasiliansk Industri / braziliai divízió / Brazilian Industry / Brasiliansk industri - SAC 0800 725 6363
Registro Anvisa: 10344420082 / 10344420371 / 10344420376 / 10344420383 / 10344420384

vide rótulo / see label / consulte la etiqueta / vedere l’etichetta / siehe Etikett / voir l’étiquette / se etiket / lásd a címkét / A se vedea eticheta / se märkning

Etkon GmbH, Lochhamer Schlag 6, 82166 Gräfelfing, Germany.

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This device is intended for a specialized procedure, which must be performed only by qualified professionals in Dental Implants for use along with the Burn-Out Coping Titanium Base or also with specific training in CAD/CAM technology. It is recommended to use the Titanium Base library on Dental Wings, exocad and 3Shape systems or softwares. For manufacturing, Neodent Digital milling centers and lab milling machines are recommended. For better results, use the product only with knowledge of the appropriate techniques and always apply them under appropriate conditions.

DESCRIPTION

The Titanium Base is an intermediary prosthetic abutment manufactured in titanium alloy according to ASTM standard F136, to be installed between implant and prosthesis.

The Titanium Base has a cylindrical format with variances of: diameter, cementable portion height and gingival height (according to the following table). It has different prosthetic interfaces in rotational (bridges) and antirotational (crown or coping) versions. It is supplied straight/without angulation and must be used in this manner.

One of its ends presents a prosthetic interface for adaptation in implants in GM, CM, HE (prosthetic platforms 3.3, 4.1/4.3 and 5.0), and II Plus prosthetic interfaces. Titanium Bases are supplied along with a screw, which can be coupled or removable. The screw from Titanium Bases AS presents a lilac coloration. The screw from the GM Titanium Base with Neo Removable Screw presents pink coloration.

Antirrotational Titanium Base (Crown or Coping): It has external retentions and three indexing guides for fitting the Crown/Coping. It is available for implants in GM, GM Exact, CM, HE, and II Plus prosthetic interfaces.

Titanium Base (for bridge): It has external retentions and is available for implants in GM, HE and II Plus prosthetic interfaces.

Prosthetic Interface	Screw	Format	Restoration Seating Diameter (mm)	Cementable Area Height (mm)	Gingival Height (mm)
CM	Coupled	Antirrotational	3.5	4.0	0.8; 1.5; 2.5; 3.5 and 4.5
HE 3.3 and II Plus	Removable	Rotational Antirrotational	4.0	4.0	0.5; 2; 3 and 4
HE 4.1/4.3			4.5	4.0	0.5; 2; 3 and 4
HE 5.0			5.3	4.0	0.5; 2; 3 and 4
GM	Removable	Rotational	3.5; 4.5 and 5.5	4.0 and 4.5	0.8; 1.5; 2.5; 3.5 and 4.5
GM AS	Removable	Antirrotational	4.0, 4.5, and 5.5	4.0 and 6.0 ²	0.8; 1.5; and 2.5
GM Exact ¹	Coupled or Removable		3.5; 4.5, 5.5, and 6.5	4.0 and 6.0 ²	0.8; 1.5; 2.5; 3.5; 4.5; and 5.5 ³
GM (Neo Removable Screw)	Removable	Rotational	3.5	4.0	0.8; 1.5; 2.5; 3.5; 4.5
			4.5 and 5.5	4.5	0.8; 1.5; 2.5; 3.5; 4.5
		Antirrotational	3.5; 4.5 and 5.5	4.0 and 6.0 ²	0.8; 1.5; 2.5; 3.5; 4.5
			6.5	4.0 and 6.0 ²	1.5; 2.5; 3.5; 4.5

¹ The GM Titanium Base with a coupled screw has a self-removal feature. Whenever counter torque is applied at the screw, the Titanium Base body simultaneously detaches from the implant. The Titanium Base of 6.5 mm diameter is supplied along with a removable screw.

² For this cementable area height, customization to 4.0 mm is allowed.

³ The GM Titanium Base with a removable screw is not available in the 5.5 mm gingival height.

INDICATIONS FOR USE

Titanium Base Abutment is a titanium abutment to be placed onto Neodent dental implants, providing support for customized prosthetic restorations such as copings, crowns or bridges.

It is indicated for cemented or screw-mounted single- or multiple-tooth cemented applications on implants installed on the maxilla or mandible.

APPLICATIONS

The Titanium Base is indicated according to the interocclusal space available, existing gingival height and three-dimensional position of the implant. The GM Exact Titanium Base 5.5 with a 0.8 mm gingival height seats on the platform and is indicated for situations in which the implant platform is at gingival height. The GM Exact Titanium Base of equal diameter to or bigger than 5.5 mm with any gingival height is applicable only for implants with diameters equal to or bigger than 5.0 mm. The GM Titanium Base with 6 mm cementable portion height allows its customization to 4 mm. The Titanium Base AS is used in single-unit prosthetic structures and presents angled access to avoid screw appearance in esthetic regions or in occlusion. They must be used in cemented or screw-retained restorations on Titanium Base which are screwed onto the implant. The HE, II Plus, and GM interfaces are suitable for single-unit or multi-unit restorations meanwhile the CM and GM Exact interfaces are indicated only for single-unit restorations. Check compatibility between the prosthetic interfaces chosen. Titanium Bases are available in GM, GM Exact, CM, HE, and II Plus prosthetic interfaces, being compatible with corresponding implants of the same lines. It is recommended to use the Titanium Base library on Dental Wings, exocad and 3Shape systems or softwares. For manufacturing, Neodent Digital milling centers and lab milling machines are recommended. In order for it to be used along with the Burn-Out Coping Titanium Base, the prosthetic structures are obtained through a conventional flow by waxing and casting.

This product's target population is presented as patients that need a rehabilitating treatment with dental implants and do not present the mentioned contraindications.

WARNING

The Titanium Bases in prosthetic interfaces GM, GM Exact, CM, HE, and II Plus (with cementable area height of 4 mm) are supplied without angulation and cannot be modified. Titanium Bases in GM and GM Exact prosthetic interfaces with cementable area height of 6 mm are supplied without angulation, but allows their customization by reducing their cementable area height to 4 mm.

The GM Rotational Titanium Base cannot be used in cases which the angulation between implants is superior than 15°.

Implants in combination with Titanium Bases with angled prosthetic structures are recommended in accordance with the following table:

Ø Implant (mm)	Mouth Region for Installation	Maximum Angulation of the Prosthetic Structure	IPS e.max CAD	3M ESPE Lava [®] plus Zirconia	CoCr ²	NIce ¹	3M Lava Esthetic ¹	Zerion LT ¹	Zerion GI ¹	Zerion UTML ¹	Zerion HTML ¹	Polycon ae ¹
3.3	1-2	17°	✓	✓	✓			✓				
	Maxilla 1	30°			✓							
	Maxilla 2 Mandible 1-2	30°	✓	✓	✓			✓				
3.5	1-2	30°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
3.75	1-5	30°	✓	✓	✓	✓	✓	✓	✓	✓		✓
≥4.0	1-8	30°	✓	✓	✓			✓	✓			✓

¹ For Brazil, these raw materials are available only for GM Titanium Bases.

² Used only for copings/ structures.

NOTE: Rotational Titanium Bases and Titanium Bases for Angled Solution (GM) can be used only with the following materials: Titanium, CoCr, Zerion LT, Zerion HTML, Polycon, and 3M Lava Plus.

HE implants with diameter 3.3 mm associated with crowns/copings/ structures angled up to 17° are recommended for incisor regions (1-2).

HE implants with diameter 3.3 mm associated with copings/structures in CoCr angled up to 30° are recommended for upper central incisor regions (Maxilla1). HE implants with diameter 3.3 mm associated with crowns / copings/ structures angled up to 30° are recommended for upper side and lower incisor regions (Maxilla 2, Mandibula 1-2).

Implants with diameter 3.5 mm associated with crowns / copings/ structures angled up to 30° are only recommended for incisor region (1-2).

Implants with diameter 3.75 mm, associated with crowns / copings/ structures angled up to 30° are only recommended for incisor, canine and premolar regions (1-5).

Implants with diameters greater than or equal to 4.0 mm associated with crowns / copings/ structures angled up to 30° are recommended for all the regions of the mouth (1-8).

CONTRAINDICATIONS

This product is contraindicated for patients exhibiting signs of allergy or hypersensitivity to the chemical ingredients of the material: Titanium alloy Ti6Al4V-ELI.

The GM Titanium Base with equal diameter to or bigger than 5.5 mm with any gingival height cannot be used with implants of smaller diameter than 5.0 mm.

HANDLING

CUSTOMIZATION OF THE CEMENTABLE AREA: The GM Exact/GM Neo Removable Screw Titanium Bases with cementable area height of 6 mm allows customization. Perform the customization in a laboratory environment with the aid of a cutting disk. Make the cut in the upper end canal. Remove all the canal area and the portion above the canal. In case of customization, maintain the minimum height of 4 mm of cementable area.

COPING/CROWN/BRIDGE PRODUCTION BY CAD/CAM TECHNOLOGY: Prepare the Coping or crown using the CAD software in accordance with the dimensions of the Titanium Base selected and the software manufacturer's instructions. Always finish the Coping/Crown before making the connection in the Titanium Base. When you use the digital platform, select the Neodent Implant library - Titanium Base with an appropriate software platform, thus facilitating the interface drawing between the Coping/Crown/Bridge and Titanium Base. This library consists of a set of 3D files containing the milling array required for the Coping/Crown/Bridge geometry.

Note: The maximum angle of the structure is variable (see table in Warning). The minimum thickness of the structure is variable depending on the material, as described in the table below:

The taper of the structure cannot exceed 8°. In the case of angled structures, the maximum overall height of the restoration cemented to the Titanium Base cannot exceed 10 mm.

Prosthetic Interface	Material	Minimum Thickness (mm)
HE / II Plus / CM / GM	CoCr	0.3
	IPS e.max CAD	0.9
	3M ESPE Lava Plus Zirconia	0.7
	Zerion LT*	0.5
CM / GM	NIce*	1.0
	3M Lava Esthetic*	0.8
	Zerion GI*	0.4
	Zerion UTML*	0.8 for anterior region 1.0 for pre-molars
	Zerion ML*	0.4 for anterior region 0.5 for posterior region
	Polycon ae*	0.6

*For Brazil, these raw materials are available only for GM Titanium Bases.

CEMENTATION: The Titanium Base must be used, in laboratory procedures, as follows: Place the Base abutment over the plaster model and pass the Fixation Screw through with a soft tightening using the Manual Screwdriver indicated (see chart).

Protect access to the screw throughout the cementation process. The cement handling must follow the manufacturer instructions. Apply the cement in the outer portion of the Titanium Base and press the restoration, fitting it according to the three available indexer guides. Press the restoration over the Titanium Base and immediately remove the excess cement. After cement setting, unscrew the infrastructure from the Analog and remove the remaining cement of the Titanium Base edge. NOTE: It is recommended to use chemically activated resin cements for bonding on metals (e.g. Panavia - Kuraray). On

Lithium Dissilicate, the use of IVOCLAR Multilink cement is required.
After these procedures, it must be handled using appropriate Drivers/Connections and adequate torque, as presented in the table below:

Prosthetic Screw model supplied along with the Titanium Base	Insertion Torque (N.cm)	Driver / Connection
Hex. Screw (HE and II Plus)	HE: 32 II Plus: 20	Hex 1.2 (mm)
CM Screw	15	Hex 0.9 (mm)
GM Screw / Neo Removable	20	Neo Short* / Neo Medium**
GM Screw for Angled Solution (AS)***	20	AS

*The Neo Manual Screwdriver (Short) is recommended for use in the laboratory phase.

**If the prosthetic restoration is cemented on the Titanium Base, it is indicated only the use of Neo Manual Screwdriver (Medium) or Neo Screwdriver Torque Connection Medium.

***After cementation, the Titanium Base AS screw gets coupled to the prosthetic abutment. Always insert the screw on Titanium Base AS before cementing the restoration.

For more information about the preparation of the restoration, please check:

- In Brazil: instructions available on www.neodent.com.br/neodent-digital;
- other countries: consult the preparation guidance from the specialized laboratory.

NOTES:

In order to use the Titanium Base AS, verify the need to install the screw in an angle using the specific Driver/Connection mentioned in the previous table.

In order for it to be used along with the Burn-Out Coping Titanium Base, follow its own instructions for use.

In order to remove the GM Titanium Base with a coupled screw, which presents a self removal feature, apply countertorque in the screw to detach it from the implant.

TRACEABILITY LABEL

This product is supplied along with three labels that allow its traceability and must be attached to the following documents: • Medical record; • Collection tax document; • Document to be delivered to the patient (ask your advisor). The identification and traceability are performed through numeric codes REF and LOT.

SANITATION

Before installation in mouth, sanitize and sterilize the final restoration. For Titanium Bases that present a removable screw, disassemble it before sanitization and sterilization.

1st step: Fully immerse the part in enzymatic detergent (diluted according to the manufacturer)
2nd step: Wash in ultrasonic cleaner for approximately 10 to 15 minutes. 3rd step: Rinse with distilled water in abundance, until completely remove the wastes from the solution. Nylon brushes use is recommended. 4th step: Dry with a clean and dry cloth or with compressed air. 5th step: Perform a visual inspection, observing if there are failures in the cleaning process. If there are still wastes, the part should be once

again immersed in detergent - 1st step - and, if necessary, the cleaning should be performed with the aid of a nylon brush. Repeat the rinsing and drying sequence.

PRESENTATION AND STERILIZATION

This product is for single only and is supplied sterile by the ethylene oxide method, being unitarily packaged. The final prosthesis must be sterilized before installation in mouth. The following sterilization methods are recommended: moist heat (steam) autoclave, gravity-displacement or dynamic-air-removal (fractionated vacuum) cycle, unwrapped, 3 minute exposure at 132 °C (270 °F). The product must be unwrapped on an appropriate tray. Use the sterilized restoration immediately after sterilization, do not store.

PRECAUTIONS

- When using the Titanium Base GM for Angled Solution with GM Scanbodies, be aware to perform scanning with the Scanbody angled/ concave side facing vestibular.
- Do not use the product if its packaging has been violated.
- Sterilization is only guaranteed if the sterile barrier (full blister) has not been damaged.
- Do not use the product beyond its expiration date.
- Inform the laboratory that any modification or damage to the prosthetic interface requires the piece to be discarded and remade.
- Ensure use of abutments that have the original Neodent prosthetic interface, which cannot be machined in the laboratory.
- The prosthetic interface of Titanium Base does not allow customization.
- Excessive wear during customization of the cementable area may compromise mechanical properties and/or invalidate the product (applicable only to Titanium Bases with cementable area height of 6 mm).
- Make sure that the manufacturer's recommendations are shared between the dentist and the laboratory during the fabrication process of the prosthesis.
- For Titanium Bases that present a removable screw, disassemble it before sanitization and sterilization.
- Regarding the systemic appearance, consider the general state of the patient's health according to the applicable literature. Regarding the local appearance, observe the conditions of the intra-oral tissues.
- This product is for single use only.
- Reuse of this product may cause: • adverse biological effects arising from product residue, microorganisms, and/or substances originating from previous uses and/or reprocessing; • changes to the original macro- and micro-structural physical, mechanical, and chemical characteristics of the product, which may affect their intended functionality. Reuse of this product does not guarantee its safety and effectiveness, and invalidates any warranty for related products.
- Inadequate surgical and/or prosthetic planning may compromise the performance of the implant/prosthesis assembly, resulting in failures of the system, such as loss or fracture of the implant, loosening of fracturing of the components and/or prosthetic screws.
- The material selection of the prosthesis

structure must consider general aspects of the patient.

- Consult the torque to be applied to the prosthetic abutment used. Excess or insufficient torque may cause undesirable results.
- For immediate load application, check the torque indication of the implant installed.
- Before each procedure, ensure a perfect fit between the pieces.
- Before each procedure, check the conditions of the Neodent instruments, always respecting their useful lifespan. Replace the instruments in case of damage, deleted markings, compromised sharpening, deformation or wear.
- Ensure that the pieces are not swallowed or inhaled by the patient.
- When choosing the Titanium Base, make sure that it is compatible with the diameter of the implant interface installed.
- During installation, make sure to align it with the insertion axis of the implant. Make sure it is perfectly seated on the Implant. To do so, it is recommended to use periapical X-rays with the parallel technique.
- Make sure to use a prosthetic screw that is compatible with the prosthetic interface and with the prosthetic abutment.
- The GM Titanium Base with a removable screw does not present a self removal feature.
- Use only the screw that is supplied along with the Titanium Base. The use of another screw can cause damage to the abutment or even to the implant, invalidating its use.
- Check passivity and perform the occlusal and interproximal adjustment after the installation of the prosthesis, taking care to avoid compromising the implant/prosthesis assembly.
- Special care must be taken during design of the prosthesis so that the interocclusal space will not be insufficient or excessive.
- Use technique and material appropriate for the molding of the implant.
- Make sure that the transfer used is compatible with the analog.
- Observe the correct selection of the analog according to the implant.
- Always use the sequence of Neodent products. The use of surgical and prosthetic instruments and or prosthetic abutments from other manufacturer does not guarantee the function of the Neodent Implant System and is exempt from the warranty of the related products.
- It is the responsibility of the dental surgeon to use the Neodent products in accordance with the instructions for use.

ADVERSE EFFECTS

Adverse effects are not expected as long as the product is used according to its instructions for use.

MAGNETIC RESONANCE IMAGING (MRI) - SAFETY INFORMATION

The Neodent Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Neodent Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need for a professional medical monitoring after the procedure and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. The professional in charge is responsible for providing these guidelines.

STORAGE CONDITIONS

This product must be stored, in its original packaging, in a clean and dry location, at a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

All the products and consumables used in the surgery for the installation of dental implants may endanger the health of anyone who handles them. Before discarding them into the environment, it is recommended to observe the current legislation and adhere to it. In case there is no current legislation, pack them in sharps collector and dispose them in medical waste.

DATE OF EXPIRATION

Written on the label.

SYMBOLS

The table describes the symbols that may be printed on the product labeling. Please, refer to the physical labeling to see which symbols are applicable to the product.

Symbols	Description
	Product size
	Catalog number
	Batch code
	Serial number
	Product manufacturing material
	UDI code
	Medical device
	Date of manufacture
	Manufacturer
	Use-by date
	Quantity
	Consult instructions for use
	Reprocessing forbidden
	Sterilized using ethylene oxide
	Sterilized using irradiation
	Do not reuse
	Do not use if package is damaged
	Do not resterilize

Symbols	Description
	Fragile
	Upper limit of temperature
	Keep dry
	Keep away from sunlight
	Authorized representative in the European Community
	Importer in the European Community
	Swiss Authorised Representative
	Single sterile barrier system
	Single sterile barrier system with protective packaging inside
	Notification required by FDA for United States market
	CE Mark
	Non-sterile

© 2021 - JJGC Indústria e Comércio de Materiais Dentários S.A. All rights reserved.
Neodent is a registered trademark of JJGC Indústria e Comércio de Materiais Dentários S.A.
Dental Wings is a trademark or registered trademark of Dental Wings Inc.
exocad is a trademark or registered trademark of exocad GmbH.
3Shape is a trademark or registered trademark of 3Shape A/S.
3M, ESPE and Lava are trademarks or registered trademarks of 3M or 3M Deutschland GmbH.
Zerion is a trademark or a registered trademark of Institut Straumann AG.
IPS e.max and Multilink are trademarks or registered trademarks of Ivoclar Vivadent AG.
Kuraray and Panavia are trademarks or registered trademarks of Kuraray Co, LTD.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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This device is intended for a specialized procedure, which must be performed only by professionals qualified in Dental Implants for use along with the Burn-Out Coping Titanium Base or also with specific training on the CAD/CAM technology. It is recommended to use the Titanium Base library in CARES Visual, Dental Wings and 3shape software. For better results, use the product knowing the appropriate techniques and always apply them under appropriate conditions.

DESCRIPTION

The Titanium Base is an intermediary prosthetic component to be installed between implant and prosthesis, being manufactured in titanium alloy according to ASTM F136 standard.

The Titanium Base has a cylindrical format with variances of: diameter, height of cementable area and gingival height (according to the following table). It has different prosthetic interfaces in rotational (bridges) and antirotational (crown or coping) versions. It is provided straight/without angulation and must be used in this manner.

One of its ends presents a prosthetic interface for adaptation on CM and GM implants. The Titanium Base is supplied with a screw, which can be coupled or removable. The screw from the AS Titanium Base presents lilac coloration. The screw from the GM Titanium Base with Neo Removable Screw presents pink coloration.

Antirotational Titanium Base (crown or coping): It has external retentions for fitting the multi-unit restoration. Available in GM, GM Exact and CM interfaces.

Rotational Titanium Base (multi-unit prosthesis): It has external retentions and three indexing guides for fitting the Crown/Coping. Available in GM interface.

Prosthetic Interface	Screw	Format	Restoration Seating Diameter (mm)	Cementable Area Height (mm)	Gingival Height (mm)
CM	Coupled	Antirotational	3.5	4.0	0.8; 1.5; 2.5; 3.5 and 4.5
GM	Removable	Rotational	3.5, 4.5 and 5.5	4.0 and 4.5	0.8; 1.5; 2.5; 3.5; 4.5 and 5.5 ³
GM AS	Removable	Antirotational	4.0, 4.5 and 5.5	4.0 and 6.0 ²	0.8, 1.5 and 2.5
GM Exact ¹	Coupled or Removable		3.5; 4.5, 5.5, and 6.5	4.0 and 6.0 ²	0.8; 1.5; 2.5; 3.5; 4.5 and 5.5 ³
GM (Neo Removable Screw)	Removable	Rotational	3.5	4.0	0.8; 1.5; 2.5; 3.5; 4.5
			4.5 and 5.5	4.5	0.8; 1.5; 2.5; 3.5; 4.5
		Antirotational	3.5; 4.5 and 5.5	4.0 and 6.0 ²	0.8; 1.5; 2.5; 3.5; 4.5
			6.5	4.0 and 6.0 ²	1.5; 2.5; 3.5; 4.5

¹ The GM Titanium Base which contains a coupled screw has a self-removal feature. Whenever counter torque is applied to the screw, the Titanium Base body simultaneously detaches from the implant. The GM Titanium Bases can also be supplied along with a removable screw.

² For this cementable height, customization to 4 mm is allowed.

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INDICATIONS FOR USE

Titanium Base Abutment is a titanium base placed onto Neodent dental implants to provide support for customized prosthetic restorations. It is used with a coping and crown, or crown alone, and is indicated for cement- or screw-retained single or multi-unit restorations single.

All digitally designed copings and/or crowns to be used with the Neodent Titanium Base Abutment System are intended to be sent to Straumann for manufacture at a validated milling center. The GM Titanium Base for Bridge is indicated for cement or screw-retained multi-unit restorations.

APPLICATIONS

The Titanium Base is recommended according to the interocclusal space available, existing gingival height and three-dimensional position of the implant. The GM Exact Titanium Base 5.5 with 0.8 gingival height seats on the platform, being indicated for situations where the implant platform is at gingival height. The GM Exact Titanium Base with diameter equal to or higher than 5.5 mm with any gingival height is applicable only for implants with diameter equal to or greater than 5.0 mm. The GM Titanium Base with 6 mm cementable height allows its customization to 4 mm. The Titanium Base AS is used in single-unit prosthetic structures and presents angled access to avoid screw appearance in esthetic regions or in occlusion. They must be used in cement- or screw-retained restorations on the Titanium Base which are screwed on the implant. They are indicated for single-unit and multi-unit restorations. Check the compatibility between the prosthetic interfaces chosen. Titanium Bases are available in the CM and GM interfaces, being compatible with corresponding implants of the same lines.

It is recommended to use the Titanium Base library in CARES Visual, Dental Wings and 3shape software.

In order for it to be used along with the Burn-Out Coping Titanium Base, the prosthetic structures are obtained through waxing and casting.

WARNING

The Titanium Bases cannot be modified.

The Ø3.5mm implants with angled abutments are recommended for incisors region only.

Small diameter implants and angled abutments are not recommended for the posterior region.

The GM Rotational Titanium Base cannot be used in cases of angulation between implants superior than 15°. The rotational Titanium Bases (GM) can be used only with the following materials: Titanium (Ticon), CoCr (Coron), Polycon ae, and Zirconia (Zerion LT).

The GM Titanium Bases for Angled Solution can be used only with the following materials: CoCr (Coron) and Zirconia (Zerion LT).

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The GM Titanium Base with diameter equal to or higher than 5.5 mm with any gingival height cannot be used with implants of smaller diameter than 5.0 mm.

HANDLING

CUSTOMIZATION OF THE CEMENTABLE AREA: The GM Exact Titanium Base with cementable area height of 6 mm allows customization. Perform the customization in a laboratory environment with the aid of a cutting disk. Make the cut in the upper end canal. Remove all the canal area and the portion above the canal. In case of customization, it is indicated to maintain the minimum height of 4 mm of cementable area.

RESTORATION BY CAD/CAM TECHNOLOGY: Prepare the Coping or crown using the CAD software in accordance with the dimensions of the Titanium Base selected and the software manufacturer's instructions. Always finish the Coping/Crown/Bridge before making the connection in the Titanium Base. Always

finish the Coping/Crown/Bridge before the Titanium Base fixation. When you use the digital platform, select the Neodent Implant kit library - Titanium Base with an appropriate software platform, thus facilitating the interface drawing between Coping/Crown/Bridge and Titanium Base. This library consists of a set of 3D files containing the milling array required for the Coping geometry. Note: The minimum wall thickness of the structure is variable depending on the material, as described in the table below. The maximum structure angulation cannot exceed 30°. The taper of the structure cannot exceed 8°. In the case of angled structures, the maximum overall height of the restoration cemented to the Titanium Base cannot exceed 10 mm.

Material	Minimum thickness (mm)
Cobalt-Chromium (Coron)	0.3
Titanium (Ticon) ¹	0.4
Zirconia (Zerion LT)	0.5
IPS e.max CAD	0.9
3M ESPE Lava Plus Zirconia	0.7
Polycon ae ²	1.0

¹ This material is only available for GM Rotational Titanium Base.

² This material is only available for CM Antirotational Titanium Base, GM Rotational Titanium Base and GM Exact Antirotational Titanium Base.

CEMENTATION: The Titanium Base must be used, in laboratory procedures, as follows: Place the Base component over the plaster model and position the Fixation Screw on it, applying a soft tight using the Screwdriver indicated (see chart). Protect access to the screw throughout the cementation process. The cement handling must follow its manufacturer instructions. Apply the cement in the outer portion of Titanium Base and press the fitting of the restoration according to the three available indexing guides. The restoration must be pressed over the Titanium Base and the excess cement removed immediately. After the cement setting, unscrew the infrastructure from the Analog and remove the surplus cement remaining in the Titanium Base edge. NOTE: It is recommended to use Panavia - Kuraray cement. On Lithium Dissilicate, the use of IVOCLAR Multilink cement is required. After these procedures, it must be handled using appropriate Drivers/Connections and adequate torque, as presented in the table below:

Prosthetic Screw model supplied along with the Titanium Base	Insertion Torque (N.cm)	Driver / Connection
CM Screw	15	Hex 0.9 (mm)
GM Screw / Neo Removable	20	Neo Short* / Neo Medium**
GM Screw for Angled Solution***	20	AS

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For more information on the preparation of the restoration, follow the preparation guidelines of the specialized laboratory.

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SANITATION

Before installation in mouth, the final restoration and the removable screw (if applicable) must be sanitized and sterilized.

1st step: Fully immerse the part in enzymatic detergent (diluted according to the manufacturer) 2nd step: Wash in ultrasonic cleaner for approximately 10 to 15 minutes. 3rd step: Rinse with distilled water in abundance, until completely remove the wastes from the solution. Nylon brushes use is recommended. 4th step: Dry with a clean and dry cloth or with compressed air. 5th step: Perform a visual inspection, observing if there are failures in the cleaning process. If there are still wastes, the part should be once again immersed in detergent - 1st step - and, if necessary, the cleaning should be performed with the aid of a nylon brush. Repeat the rinsing and drying sequence.

PRESENTATION AND STERILIZATION

This product is for single use and supplied sterile, being unitarily packaged.

Sterilization method: ethylene oxide. This product is unitarily provided, conditioned in blister (surgical grade film and paper type) packing. The final restoration must be sterilized before installing in mouth. The following sterilization methods are recommended: moist heat (steam) autoclave, gravity-displacement or dynamic-air-removal (fractionated vacuum) cycle, unwrapped, 3-minute exposure at 132 °C (270 °F). The product must be unwrapped on an appropriate tray. Use the sterilized restoration immediately after sterilization, do not store.

PRECAUTIONS

- When using the Titanium Base GM for Angled Solution with GM Scanbodies, be aware to perform scanning with the Scanbody angled/ concave side facing vestibular.
- Do not use the product if its packaging is damaged.
- Sterilization is only guaranteed if the packaging has not been damaged.
- Do not use the product with its validity expired.
- Inform the laboratory that any modification or damage to the prosthetic interface requires that the piece to be discarded and remade.
- Ensure that only components with the original Neodent prosthetic interface are being used, which cannot be machined in the laboratory.
- The prosthetic interface of Titanium Base does not allow customization.
- Excessive wear during the customization of the cementable area can compromise the mechanical properties and/or invalidate the product (applicable only to Titanium Bases with cementable area height of 6 mm).
- Make sure that the manufacturer's recommendations are shared between the dentist and the laboratory during the fabrication process of the prosthesis.
- Regarding the systemic aspect, consider the general state of the patient's health according to the applicable literature. Regarding the local appearance, observe the conditions of the intra-oral tissues.
- This product is for single use only.
- Reuse of this product may cause: • adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing; • changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put at risk the desired functionality. The reuse of this product does not guarantee its safety or efficacy and disclaims any warranty of products.
- Inadequate surgical and/or prosthetic planning may compromise the performance of the implant/ prosthesis assembly, resulting in failures of the system, such as loss or fracture of the implant, loosening or fracturing of abutments and/or prosthetic screws.
- The material selection of the prosthesis structure must consider general aspects of the patients.
- Consult the torque to be applied to the prosthetic abutment chosen. Excess or insufficient torque may cause undesirable results.
- For immediate load application, check the torque indication of the implant installed.
- Before each procedure, ensure a perfect fit between the pieces.

- Before each procedure, check the conditions of the Neodent instruments, always respecting their useful lifespan. Replace the instruments in case of damage, deleted markings, compromised sharpening, deformations, or wear.
- Ensure that the pieces are not swallowed or inhaled by the patient.
- When choosing the Titanium Base, make sure that it is compatible with the diameter of the implant interface installed.
- During installation, make sure to align it with the insertion axis of the implant. Make sure it is perfectly seated on the Implant. To do so, it is recommended to use periapical X-rays with the parallel technique.
- Make sure to use a prosthetic screw that is compatible with the prosthetic interface and with the prosthetic abutment.
- The GM Titanium Base with a removable screw does not present a self-removal feature.
- Use only the screw that is supplied along with the Titanium Base. The use of another screw can cause damage to the abutment or even to the implant, rendering them unusable.
- Check passivity and perform the occlusal and interproximal adjustment after the installation of the prosthesis, avoiding impairment of the implant/prosthesis assembly.
- Special care must be taken during design of the prosthesis so that the interocclusal space will not be insufficient or excessive.
- Use appropriate technique and material for the molding of the implant.
- Make sure that the transfer used is compatible with the analog.
- Observe the correct selection of the analog according to the implant.
- Always use the sequence of Neodent products. The use of surgical and prosthetic instruments and or prosthetic abutments from other manufacturer does not guarantee the function of the Neodent Implant System and exempts any warranty of related products.
- It is the dental surgeon's responsibility to use Neodent products according to their instructions for use.

ADVERSE EFFECTS

Adverse effects are not expected as long as the product is used in accordance with its instructions for use.

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This product is manufactured from a metallic material that can be affected in magnetic resonance field.

MR Conditional

Non-clinical testing and magnetic resonance (MRI) simulations have shown that Neodent implantable devices manufactured in metallic material are MR Conditional. A patient that has a Neodent implantable medical device can be scanned safely under the following conditions:

- Static magnetic field of 1.5T and 3T, only.
- Maximum spatial gradient magnetic field of 4,000 gauss/cm (40 T/m).
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg and head average SAR of 3.2 W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the

Normal Operating Mode.

Under the scan conditions defined, the implants, as well as clinically relevant implant constructs, from the Neodent Implant System are expected to produce a maximum temperature rise of 4.9°C after 15-minutes of continuous scanning (i.e., per pulse sequence). In non-clinical testing, the image artifact caused by the implants from the Neodent Implant System extend approximately 10 mm from this device when imaged with a gradient echo pulse sequence and a 3T MR system.

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Symbols	Description
	Product size
	Catalog number
	Batch code
	Serial number
	Product manufacturing material
	UDI code
	Medical device
	Date of manufacture
	Manufacturer
	Use-by date
	Quantity
	Consult instructions for use

Symbols	Description
	Reprocessing forbidden
	Sterilized using ethylene oxide
	Sterilized using irradiation
	Do not reuse
	Do not use if package is damaged
	Do not resterilize
	Fragile
	Upper limit of temperature
	Keep dry
	Keep away from sunlight
	Authorized representative in the European Community
	Importer in the European Community
	Swiss Authorised Representative
	Single sterile barrier system
	Single sterile barrier system with protective packaging inside
	Notification required by FDA for United States market
	CE Mark
	Non-sterile

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Neodent is a registered trademark of JJGC Indústria e Comércio de Materiais Dentários S.A. Dental Wings is a trademark or a registered trademark of Dental Wings Inc. 3Shape is a trademark or a registered trademark of 3Shape A/S. 3M, ESPE and Lava are trademarks or registered trademarks of 3M or 3M Deutschland GmbH. Multilink and IPS e.max are trademarks or registered trademarks of Ivoclar Vivadent AG. Kuraray and Panavia are trademarks or registered trademarks of Kuraray Co, LTD. Straumann is a trademark or registered trademark of Institut Straumann Holding AG.

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