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DESCRIPTION

The CoCr Abutment For Crown Set is composed of Base, Screw, and Analog.
The Base is composed of two parts: coping and base. The coping has two channels for wax retention during the crown confection and is manufactured in polymer (POM). The Base has 1 mm of transmucosal height, interface compatible with the GM implant and is manufactured in Cobalt-Chromium-based alloy.
The Analog presents a geometry with a channel that favors its retention in the plaster model and is compatible with the GM interface. It is available in yellow.
The Screw is formed by two different parts: head and body. The head has a geometry to fit into the Neo Screwdriver / Neo Connection that will enable its use and seating in the internal part of the Base. The body has a thread with 1.6 mm diameter.
The Screw and the Analog from the CoCr Abutment for Crown Set are of laboratory use and manufactured in titanium alloy.
For definitive use, use the Neotorque Screw.

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APPLICATIONS

The CoCr Abutment For Crown Set is composed of Base, Screw, and Analog.
The Base is indicated for the confection of unitary definitive prostheses that are screwed or cemented on the GM Implant. The coping is used in the laboratory phase for fabricating the prosthesis on plaster models and implant analog. Upon foundry, it becomes part of the definitive prosthesis. The Base has compatible interface with the GM implant and cannot be modified.
The Analog is to be used in the laboratory phase while obtaining the plaster model for fabricating the prosthesis on implant, reproducing the dimensions and positioning of the implant.
The Screw is to be used in the laboratory phase for the fixation of the Base on the Analog, fixed in the plaster model.

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CONTRAINDICATIONS

The CoCr Abutment for Crown Set is contraindicated in multiple unitary rehabilitations and on implants placed subcrestally.
This product is provided straight/without angulation and must be used in this format.
This product is contraindicated for patients exhibiting signs of allergy or hypersensitivity to the chemical elements of the material: Cobalt-Chromium-based alloy.
This product is contraindicated for insufficient interocclusal space and unsatisfactory three-dimensional position of the implant.

HANDLING

For use in two-stage procedures, a prior preparation of the soft tissues can be performed with the use of a Healing Abutment specific for the CoCr Abutment for Crown Set.
Perform the molding of the implant by using the Impression Coping for GM implant according to proper techniques.

LABORATORY PHASE

After obtaining the mold, screw the Analog into the Impression Coping used. Apply maximum torque of 20 N.cm. Pour rubbery material to simulate the soft tissue of the patient. Pour the plaster on the mold, covering the whole Analog. After the plaster setting phase, separate the mold from it. Unscrew/disassemble the Impression Coping from the Analog.
Screw the Base on the Analog by using the Neo Connection / Neo Screwdriver, with the recommended torque of 20 N.cm. Customize the coping as needed. For customization, the minimum interocclusal height of 5 mm must be maintained. Perform the closing.
The obtained standard must be included in a compatible coating with the alloy that will be used in the casting process. The chosen alloy must be Cobalt-Chromium-based and compatible with the aesthetic material that is planned for the prosthesis. It is recommended to use phosphate coating to include the standard in wax. To obtain infrastructures for metal ceramics prosthesis, a compatible alloy must be chosen. It is important that the ceramics cannot be applied directly on the Base and, therefore, the waxing must ensure coverage of the entire external area of the abutment on which ceramics will be applied.
The application of ceramics (specific for alloys of this kind) directly over an area uncovered by overcasting metal alloy may cause cracks. Special care must be taken during waxing at the edges of the abutment to avoid invasion at the edges that will come into contact with the implant. The thickness of

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CLINICAL PHASE

Use the Neotorque Screw for the intra oral fixation of the Base after casting. Passivity and adaptation tests of the prosthesis structure must be performed.
Note: the Screw (within this set) is of laboratory use. The Neotorque Screw for intra oral fixation must be bought separately.

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SANITIZATION

Before installation in mouth, sanitize and sterilize the final abutment.
1st step: Fully immerse the part in enzymatic detergent (diluted according to the manufacturer)
2nd step: Wash in ultrasonic cleaner for approximately 10 to 15 minutes. 3rd step: Rinse with distilled water in abundance, until completely remove the wastes from the solution. Nylon brushes use is recommended. 4th step: Dry with a clean and dry cloth or with compressed air. 5th step: Perform a visual inspection, observing if there are failures in the cleaning process. If there are still wastes, the part should be once again immersed in detergent – 1st step – and, if necessary, the cleaning should be performed with the aid of a nylon brush. Repeat the rinsing and drying sequence.

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PRECAUTIONS

- Be sure to use the corresponding set with the planned implant according to use indication. Its parts cannot be used separately.
- Cleaning and sterilization of the final abutment before its installation in mouth is recommended.
- Before installing the product, make sure it has the same prosthetic interface as the implant. Ensure that the stability of the implant is sufficient to withstand the installation torque of the prosthetic abutment and functional load, according to the instructions for use of the implant.
- WARNING: For overcasting, always use

Cobalt-Chromium-based alloy with lower melting temperature. Other metal alloys do not guarantee proper bonding. Ceramics cannot be applied directly onto the prosthetic abutment fitting strap. Melting must be made in the range from 1120 to 1150°C (2048 – 2102°F).

- Inadequate surgical and/or prosthetic planning can compromise the performance of the implant / prosthesis assembly, resulting in system failure, such as loss or fracture of the implant, loosening or fracture of components and/or prosthetic screws.
- The material selection of the prosthesis structure must consider general aspects of the patient.
- This product is for single use only and cannot be reprocessed.
- Reuse of this product may cause:
 - adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing;
 - changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put at risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and exempts any warranty of products.
- Do not use the product if the packaging is damaged.
- Do not use the product with its validity expired.
- For immediate load application, check the torque indication of the implant installed.
- See the torque to be given on the prosthetic abutment chosen. The excess or lack of torque can lead to undesirable results.
- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.
- During installation, make sure to align it with the insertion axis of the implant. Make sure it is perfectly seated on the Implant. To do so, it is recommended to use periapical X-rays with the parallel technique.
- Make sure to use a prosthetic screw that is compatible with the prosthetic interface and with the prosthetic abutment.
- Check passivity and perform occlusal and interproximal adjustment after installation of the prosthesis, avoiding impairment of the implant/ prosthesis assembly.
- Before each procedure check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is damage, deleted markings, compromised sharpening, deformation or wear.
- Always use the Neodent product sequence. The use of prosthetic abutments and/or instruments from other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty.
- It is the dental surgeon's responsibility to use the Neodent products according to the instructions for use.

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Symbols	Description
	Swiss Authorised Representative
	Single sterile barrier system
	Single sterile barrier system with protective packaging inside
	Notification required by FDA for United States market
	CE Mark
	Non-sterile
	Contains hazardous substances

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Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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This device is intended for a specialized procedure, which must be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques and always apply them under appropriate conditions, including in an operating room.

DESCRIPTION

The product package contains the GM Exact CoCr Abutment for Crown, abutment Screw, and implant Analog.

The GM Exact CoCr Abutment for Crown is composed of two parts: the CoCr alloy that interfaces directly to the implant and a polymer cylinder part for fabrication of a cast prosthesis by a burn-out technique. The cylinder part has two channels for wax retention during the crown fabrication and is manufactured in polyoxymethylene polymer (POM). The CoCr alloy component has 1 mm of transmucosal height, and the interface is compatible with the GM implant.

The implant Analog has a geometry with a channel for retention in the plaster model and is compatible with the GM interface. It is available in a yellow color.

The abutment Screw is formed by two different parts: head and body. The head has a geometry to fit into the Neo Screwdriver / Neo Connection that will enable its use and seating in the internal part of the GM Exact CoCr Abutment for Crown. The Screw body has a thread with 1.6 mm diameter.

The abutment Screw must be cleaned and sterilized before its use, being used as a definitive screw afterwards.

The implant Analog is for laboratory use only and is manufactured in titanium alloy. For laboratory use, the working screw must be bought separately.

INDICATIONS FOR USE

The Neodent Implant System is intended to be surgically placed in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, to restore chewing function. It may be used with single-stage or two-stage procedures, for single or multiple unit restorations, and may be loaded immediately when good primary stability is achieved and with appropriate occlusal loading.

APPLICATIONS

The GM Exact CoCr Abutment for Crown is indicated for the fabrication of single-unit definitive prosthesis that is screwed or cemented on the GM Implant. The GM Exact CoCr Abutment for Crown is used in the laboratory phase for fabricating the prosthesis on plaster models and implant analog. Upon foundry, it becomes part of the definitive prosthesis. The CoCr alloy part has an interface compatible with the GM implant and the interface cannot be modified.

The implant Analog is indicated for use in the laboratory phase while obtaining the plaster model for fabricating the prosthesis on implant, reproducing the dimensions and position of the implant.

After being cleaned and sterilized, the abutment Screw is indicated for definitive use.

Prosthetic Interface	Item	Implant Diameter (mm)	Gingival Height (mm)
GM Exact	CoCr Abutment For Crown Set	3.5 and 3.75	1
		4.0 and 4.3	
		5.0 and 6.0	

Check compatibility between the prosthetic interfaces chosen. The GM Exact CoCr Abutment for Crown is available in the GM prosthetic interface and is exclusively compatible with GM implants.

CONTRAINDICATIONS

The GM Exact CoCr Abutment for Crown is contraindicated in multiple unit rehabilitations and on implants placed subcrestally.

This product is provided straight/without angulation and must be used in this manner.

This product is contraindicated for patients exhibiting signs of allergy or hypersensitivity to the chemical elements of the material: Cobalt-Chromium-based alloy.

This product is contraindicated for insufficient interocclusal space and unsatisfactory three-dimensional position of the implant.

HANDLING

For use in two-stage procedures, a prior preparation of the soft tissues can be done with the use of a Healing Cap specific for GM line.

Perform the molding of the implant by using the Impression Coping for GM implant according to appropriate techniques.

LABORATORY PHASE

After obtaining the mold, screw the Analog into the Impression Coping. Pour rubbery material to simulate the soft tissue of the patient. Pour the plaster on the mold, covering the whole Analog. After the plaster setting phase, unscrew/unplug the Impression Coping from the Analog and remove (separate) the mold from the plaster model.

Screw the GM Exact CoCr Abutment for Crown on the Analog with the Neo Connection / Neo Screwdriver, to the recommended torque of 20 N.cm. Prep the POM cylinder part of the GM Exact CoCr Abutment for Crown as needed. For prepping, the minimum interocclusal height of 5 mm must be maintained. Perform waxing.

The wax pattern obtained over GM Exact CoCr Abutment for Crown must be submerged in a coating material compatible with the alloy that will be used in the casting process. The chosen alloy must be Cobalt-Chromium-based and compatible with the aesthetic material that is planned for the prosthesis. For overcasting, always use Cobalt-Chromium-based alloy with a low melting temperature. Other metal

alloys do not guarantee proper bonding. Ceramics cannot be applied directly onto the prosthetic abutment fitting strap. Melting must be made in the range from 1120 – 1150 °C (2048 – 2102 °F).

The use of a phosphated coating is recommended as a material to submerge the wax pattern.

To obtain infrastructures for metal ceramics prosthesis, a compatible alloy must be chosen. It is important that the ceramics cannot be applied directly on the GM Exact CoCr Abutment for Crown and, therefore, the waxing must ensure a coating over the entire external area of the abutment on which ceramics will be applied.

The ceramic application (specific for alloys of this type) directly over an uncovered area by metal alloy overcasting may cause cracks. Special care must be taken during waxing at the edges of the abutment to avoid invasion at the edges that will come into contact with the implant.

The thickness of the waxing must be, at least, of 0.5 mm and may be reduced to 0.3 mm after overcasting. After casting, observe the area corresponding to the prosthetic screw seat.

The original shape must be preserved without the presence of positive metal bubbles. During finishing and polishing procedures, the prosthetic interface cannot be damaged. During the alloy finishing process, do not use corrosive materials, such as those containing iron particles. It is also recommended to machine the alloy at low speed. Do not modify in any way the CoCr alloy abutment that interfaces directly to the implant.

The CoCr Abutment for Crown is not intended to be placed or cast/modified to provide an angle/divergence correction.

- GM Exact CoCr Abutment for Crown is not intended to be placed or casted/modified to provide an angle/divergence correction.

- GM Exact CoCr Abutment for Crown is for single-unit use only.

- Casting material = cobalt-chromium based alloy compatible with the aesthetic material that is planned for the prosthesis (for example ISO 5832-4 material).

- For cementation, Panavia™ – Kuraray cement is recommended as luting material.

Design Parameters for the final cast abutment:

Parameters	Minimum	Maximum
Diameter	4.1 mm	11.5 mm
Wall thickness	0.3 mm	4.1 mm
Post height	4.0 mm	12.0 mm
Gingival height	1.0 mm	1.0 mm
Taper of cast-to portion	0°	4°

CLINICAL PHASE

Use the abutment Screw for the intra oral fixation of the final prosthesis. Passivity and adaptation tests of the prosthesis structure must be performed.

Note: the abutment Screw within the product package set is of definitive use after being cleaned and sterilized, as described in the sequence. The working screw for laboratory use must be acquired separately.

TRACEABILITY LABEL

This product is supplied along with three labels that allow its traceability and must be attached to the following documents: • Medical record; • Collection tax document; • Document to be

delivered to the patient (ask your advisor). The identification and traceability are performed through numeric codes REF and LOT.

SANITATION

Before installation in mouth, you must sanitize and sterilize the final restoration and the abutment screw.

1st step: Fully immerse the part in enzymatic detergent (diluted according to the manufacturer)
2nd step: Wash in ultrasonic cleaner for approximately 10 to 15 minutes. 3rd step: Rinse with distilled water in abundance, until completely remove the wastes from the solution. Nylon brushes use is recommended. 4th step: Dry with a clean and dry cloth or with compressed air. 5th step: Perform a visual inspection, observing if there are failures in the cleaning process. If there are still wastes, the part should be once again immersed in detergent – 1st step – and, if necessary, the cleaning should be performed with the aid of a nylon brush. Repeat the rinsing and drying sequence.

PRESENTATION AND STERILIZATION

This product package set is for single use and supplied non-sterile, being unitarily packaged. The final prosthetic component and abutment Screw must be sterilized before installing in mouth.

The following sterilization methods are recommended: moist heat (steam) autoclave, gravity-displacement or dynamic-air-removal (fractionated vacuum) cycle, unwrapped, 3-minute exposure at 132 °C (270 °F). The product must be unwrapped on an appropriate tray. Use the sterilized restoration immediately after sterilization, do not store.

Caution:

It is the responsibility of the user to establish whether their sterilizer has been cleared by the FDA to meet these recommended parameters, and to use accessories (BIs, CIs, and wraps/pouches/containers) cleared by FDA and labeled for use with the recommended sterilization parameters.

PRECAUTIONS

- Be sure to use the corresponding set with the planned implant as indicated. Assembly parts cannot be used separately.
- Before installing the product, make sure it has the same prosthetic interface as the implant. Ensure that the stability of the implant is sufficient to support the installation torque of the prosthetic abutment and functional load, according to the directions for use of the implant.
- **WARNING:** For overcasting, always use Cobalt-Chromium-based alloy with a lower melting temperature. Other metal alloys do not guarantee proper bonding. Ceramics cannot be applied directly onto the prosthetic abutment fitting strap. Melting must be made in the range from 1120 to 1150 °C (2048 – 2102 °F).
- Inadequate surgical and/or prosthetic planning can compromise the performance of the implant / prosthesis assembly, resulting in system failure, such as loss or fracture of the implant, loosening or fracture of components and/or prosthetic screws.
- The material selection of the prosthesis structure must consider general aspects of the patient.
- This product is for single use only and cannot be reprocessed.
- Reuse of this product may cause:
 - adverse

biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing; • changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put at risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.

- Do not use the product if the packaging is damaged.
- Do not use the product after its expiry date.
- For immediate load application, check the torque indication of the implant installed.
- See the torque to be given on the prosthetic abutment chosen. The excess or lack of torque can lead to undesirable results.
- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.
- During installation, make sure to align it with the insertion axis of the implant. Make sure it is perfectly seated on the Implant. To do so, it is recommended to use periapical X-rays with the parallel technique.
- Make sure to use a prosthetic screw that is compatible with the prosthetic interface and with the prosthetic abutment.
- Check passivity and perform occlusal and interproximal adjustment after installation of the prosthesis, avoiding impairment of the implant/ prosthesis assembly.
- Before each procedure, check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is any damage, deleted markings, compromised sharpening, deformation or wear.
- Always use the Neodent product sequence. The use of prosthetic abutments and/or instruments from other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty.
- It is the professional's responsibility to use the Neodent products according to the instructions for use.

ADVERSE EFFECTS

Adverse effects are not expected as long as the product is used in accordance with its instructions for use.

MAGNETIC RESONANCE IMAGING (MRI) – SAFETY INFORMATION

The Neodent Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Neodent Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need for a professional medical monitoring after the surgery and to obey the guidelines regarding precautions, hygiene and prescription of drugs. The professional in charge is responsible for providing these guidelines.

STORAGE CONDITIONS

This product must be stored, in its original packaging, in a clean and dry location, at a maximum temperature of 40 °C and protected from direct sunlight.

DISPOSAL OF MATERIAL

All the products and consumables used in the surgery for the installation of dental implants may endanger the health of anyone who handles them. Before discarding them into the environment, it is recommended to observe the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

SYMBOLS

The table describes the symbols that may be printed on the product labeling. Please, refer to the physical labeling to see which symbols are applicable to the product.

Symbols	Description
	Product size
	Catalog number
	Batch code
	Serial number
	Product manufacturing material
	UDI code
	Medical device
	Date of manufacture
	Manufacturer
	Use-by date
	Quantity
	Consult instructions for use
	Reprocessing forbidden
	Sterilized using ethylene oxide
	Sterilized using irradiation
	Do not reuse
	Do not use if package is damaged
	Do not resterilize
	Fragile
	Upper limit of temperature
	Keep dry
	Keep away from sunlight
	Authorized representative in the European Community
	Importer in the European Community

Symbols	Description
	Swiss Authorised representative
	Single sterile barrier system
	Single sterile barrier system with protective packaging inside
Rx only	Notification required by FDA for United States market
	CE Mark
NON-STERILE 	Non-sterile
	Contains hazardous substances

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This device is intended for a specialized procedure, which must be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques and always apply them under appropriate conditions.

DESCRIPTION

The CoCr Abutment For Crown Set is composed of Base, Screw, and Analog.

The Base is composed of two parts: coping and base. The coping has two channels for wax retention during the crown confection and is manufactured in polymer (POM). The Base has 1 mm of transmucosal height, interface compatible with the GM implant and is manufactured in Cobalt-Chromium-based alloy.

The Analog presents a geometry with a channel that favors its retention in the plaster model and is compatible with the GM interface. It is available in yellow.

The Screw is formed by two different parts: head and body. The head has a geometry to fit into the Neo Screwdriver / Neo Connection that will enable its use and seating in the internal part of the Base. The body has a thread with 1.6 mm diameter.

The Screw and the Analog from the CoCr Abutment for Crown Set are of laboratory use and manufactured in titanium alloy.

For definitive use, use the Neotorque Screw.

INDICATIONS FOR USE

Indicated for the confection of definitive prosthesis to obtain prosthesis infrastructure on implants.

APPLICATIONS

The CoCr Abutment For Crown Set is composed of Base, Screw, and Analog.

The Base is indicated for the confection of unitary definitive prostheses that are screwed or cemented on the GM Implant. The coping is used in the laboratory phase for fabricating the prosthesis on plaster models and implant analog. Upon foundry, it becomes part of the definitive prosthesis. The Base has compatible interface with the GM implant and cannot be modified.

The Analog is to be used in the laboratory phase while obtaining the plaster model for fabricating the prosthesis on implant, reproducing the dimensions and positioning of the implant.

The Screw is to be used in the laboratory phase for the fixation of the Base on the Analog, fixed in the plaster model.

Prosthetic Interface	Set	Implant Diameter	Gingival Height
GM Exact	CoCr Abutment For Crown Set 3.5/3.75	3.5 and 3.75 mm	1 mm
	CoCr Abutment For Crown Set 4.0/4.3	4.0 and 4.3 mm	
	CoCr Abutment For Crown Set 5.0/6.0	5.0 and 6.0 mm	

Check compatibility between the prosthetic interfaces chosen. The CoCr Abutment for Crown Set is available in the GM prosthetic interface and is exclusively compatible with Grand Morse implants.

CONTRAINDICATIONS

The CoCr Abutment for Crown Set is contraindicated in multiple unitary rehabilitations and on implants placed subcrestally.

This product is provided straight/without angulation and must be used in this format.

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HANDLING

For use in two-stage procedures, a prior preparation of the soft tissues can be performed with the use of a Healing Abutment specific for the CoCr Abutment for Crown Set.

Perform the molding of the implant by using the Impression Coping for GM implant according to proper techniques.

LABORATORY PHASE

After obtaining the mold, screw the Analog into the Impression Coping used. Apply maximum torque of 20 N.cm. Pour rubbery material to simulate the soft tissue of the patient. Pour the plaster on the mold, covering the whole Analog. After the plaster setting phase, separate the mold from it. Unscrew/disassemble the Impression Coping from the Analog.

Screw the Base on the Analog by using the Neo Connection / Neo Screwdriver, with the recommended torque of 20 N.cm. Customize the coping as needed. For customization, the minimum interocclusal height of 5 mm must be maintained. Perform the closing.

The obtained standard must be included in a compatible coating with the alloy that will be used in the casting process. The chosen alloy must be Cobalt-Chromium-based and compatible with the aesthetic material that is planned for the prosthesis. It is recommended to use phosphate coating to include the standard in wax. To obtain infrastructures for metal ceramics prosthesis, a compatible alloy must be chosen. It is important that the ceramics cannot be applied directly on the Base and, therefore, the waxing must ensure coverage of the entire external area of the abutment on which ceramics will be applied.

The application of ceramics (specific for alloys of this kind) directly over an area uncovered by overcasting metal alloy may cause cracks. Special care must be taken during waxing at the edges of the abutment to avoid invasion at the edges that will come into contact with the implant. The thickness of

the waxing must be, at least, of 0.5 mm and may be reduced to 0.3 mm after overcasting. After casting, observe the area corresponding to the prosthetic screw seat. The original shape must be preserved without the presence of positive metal bubbles. During finishing and polishing procedures, the prosthetic interface cannot be damaged. During the alloy finishing process, do not use corrosive materials, such as those containing iron particles. It is also recommended to machine the alloy at low speed.

CLINICAL PHASE

Use the Neotorque Screw for the intra oral fixation of the Base after casting. Passivity and adaptation tests of the prosthesis structure must be performed.

Note: the Screw (within this set) is of laboratory use. The Neotorque Screw for intra oral fixation must be bought separately.

TRACEABILITY LABEL

This product is supplied along with three labels that allow its traceability and must be attached to the following documents: • Medical record; • Collection tax document; • Document to be delivered to the patient (ask your advisor). The identification and traceability are performed through numeric codes REF and LOT.

SANITIZATION

Before installation in mouth, sanitize and sterilize the final abutment.

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PRESENTATION AND STERILIZATION

This product is for single use and supplied non-sterile, being unitarily packaged. The final abutment must be sterilized before its installation in mouth. The following sterilization methods are recommended: moist heat (steam) autoclave, gravity-displacement or dynamic-air-removal (fractionated vacuum) cycle, unwrapped, 3 minute exposure at 132 °C (270 °F). The product must be unwrapped on an appropriate tray. Use the sterilized restoration immediately after sterilization, do not store.

PRECAUTIONS

- Be sure to use the corresponding set with the planned implant according to use indication. Its parts cannot be used separately.
- Cleaning and sterilization of the final abutment before its installation in mouth is recommended.
- Before installing the product, make sure it has the same prosthetic interface as the implant. Ensure that the stability of the implant is sufficient to withstand the installation torque of the prosthetic abutment and functional load, according to the instructions for use of the implant.

- **WARNING:** For overcasting, always use Cobalt-Chromium-based alloy with lower melting temperature. Other metal alloys do not guarantee proper bonding. Ceramics cannot be applied directly onto the prosthetic abutment fitting strap. Melting must be made in the range from 1120 to 1150°C (2048 – 2102°F).
- Inadequate surgical and/or prosthetic planning can compromise the performance of the implant / prosthesis assembly, resulting in system failure, such as loss or fracture of the implant, loosening or fracture of components and/or prosthetic screws.
- The material selection of the prosthesis structure must consider general aspects of the patient.
- This product is for single use only and cannot be reprocessed.
- Reuse of this product may cause:
 - adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing;
 - changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put at risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and exempts any warranty of products.
- Do not use the product if the packaging is damaged.
- Do not use the product with its validity expired.
- For immediate load application, check the torque indication of the implant installed.
- See the torque to be given on the prosthetic abutment chosen. The excess or lack of torque can lead to undesirable results.
- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.
- During installation, make sure to align it with the insertion axis of the implant. Make sure it is perfectly seated on the Implant. To do so, it is recommended to use periapical X-rays with the parallel technique.
- Make sure to use a prosthetic screw that is compatible with the prosthetic interface and with the prosthetic abutment.
- Check passivity and perform occlusal and interproximal adjustment after installation of the prosthesis, avoiding impairment of the implant/ prosthesis assembly.
- Before each procedure check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is damage, deleted markings, compromised sharpening, deformation or wear.
- Always use the Neodent product sequence. The use of prosthetic abutments and/or instruments from other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty.
- It is the dental surgeon's responsibility to use the Neodent products according to the instructions for use.

SAFETY INFORMATION

Cobalt is classified as a hazardous substance. Neodent performed an assessment of risks and benefits of the cobalt present in Neodent CoCr alloy. Overall, according to the assessment, the material is safe, as long as the number of inlays does not exceed 3 units per patient. However, in the interest of your patient, consider alternative solutions.

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Instruct the patient as to the need for a professional medical monitoring after the procedure and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. The professional in charge is responsible for providing these guidelines.

STORAGE CONDITIONS

This product must be stored, in its original packaging, in a clean and dry location, at a maximum temperature of 40°C and protected from direct sunlight.

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Symbols	Description
	Product size
	Catalog number
	Batch code
	Serial number
	Product manufacturing material
	UDI code
	Medical device
	Date of manufacture
	Manufacturer
	Use-by date
	Quantity
	Consult instructions for use
	Reprocessing forbidden
	Sterilized using ethylene oxide
	Sterilized using irradiation
	Do not reuse

Symbols	Description
	Do not use if package is damaged
	Do not re sterilize
	Fragile
	Upper limit of temperature
	Keep dry
	Keep away from sunlight
	Authorized representative in the European Community
	Importer in the European Community
	Swiss Authorised Representative
	Single sterile barrier system
	Single sterile barrier system with protective packaging inside
	Notification required by FDA for United States market
	CE Mark
	Non-sterile
	Contains hazardous substances

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© 2021 — JJGC Indústria e Comércio de Materiais Dentários S.A. All rights reserved.
© 2021 – JJGC Indústria e Comércio de Materiais Dentários S.A. All rights reserved.
Panavia™ is a trademark of Kuraray Co., LTD, Japan.
© 2022 – JJGC Indústria e Comércio de Materiais Dentários S.A. All rights reserved.
© 2022 – JJGC Indústria e Comércio de Materiais Dentários S.A. Vse pravice pridržane.